

# **SPORT SAFE Testing Service, Inc.**

Founded 1996 481 Village Park Dr. & Powell, OH 43065 & (614) 847-0847 & Fax (614) 847-0874

## **Toxicology Testing vs. Oral Fluids**

## **Enzyme Immunoassay (EIA)**



### **Testing Methodology**

All drugs of abuse screens are performed using enzyme immunoassay (EIA). Confirmation of a positive drug screen are performed by gas chromatography/ mass spectrometry (GC/MS) which conclusively identifies a drug at no additional expense (from SPORT SAFE).

## **Detectable Drugs:**

Alcohol\* Marijuana\*
Amphetamines\* MDMA\* (Ecstacy)
Anabolic Steroids Methadone\*
Barbiturates\* Nicotine
Benzodiazepines\* Opiates\* (Codeine)
Cocaine\* Phencyclidine\*
LSD Propoxyphene\*

\*11 Drug Test Panel

#### **Split Specimen Possible**

**Detection Levels:** Screening Confirmation Marijuana (metabolite) 20 ng/mL 8 ng/mL

**Detection Time (THC):** Occasional Use 7-10 days Chronic Use: 30+ days

**Accuracy:** 0% False Positives. Results upheld by U.S. Supreme Court many times.

**Adulteration Testing**: 17 separate test parameters

**Cost**: 10 Drug Test Panel \$26.00 each Ecstasy \$3.00 additional Includes MRO Services

#### **Oral Fluids**

## **Testing Methodology**

Limited drugs are screened for by first inserting a pad between the cheek and gum for 2-5 minutes. The



specimen is then shipped to the lab for evaluation.

Up to 8 drugs of abuse can be screened for using enzyme immunoassay (EIA). Confirmation of a positive drug screen are performed by (GC/MS) which conclusively identifies a drug, *at additional expense*.

Oral fluid tests look for "parent drug compound for THC" and not a metabolite. Parent drug compound can be washed from the mouth very quickly after use.

#### **Detectable Drugs:**

Amphetamine Marijuana (parent compound)

Barbiturates Methamphetamine

Benzodiazepines Opiates Cocaine PCP

#### **Split Specimen NOT possible!**

Detection Levels: Screening Only

Marijuana (THC only) 1 ng/mL

**Detection Time (THC):** ONLY 10 to 24 hours

**Accuracy**: "After 48 to 72 hours oral fluids are less effective at detecting positives than urine.."

Results have not been upheld in any U.S. Court.

M. Peat, Ph.D.; LabOne, Inc.. April 12, 2000

Adulteration Testing: None

**Cost**: Intercept<sup>TM</sup> Oral Collection \$25-40.00 each Confirmation extra! MRO Services Extra!